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Mask Cleaning and Disinfection procedure

This cleaning and disinfection procedure applies to the following masks:
F300, N300, N500 series masks.

1 Precautions

- Failure to follow the manufacturer's specified cleaning agents, disinfectants, or cleaning/disinfection methods, or exceeding the maximum recommended reprocessing cycles may adversely affect the mask, ultimately compromising treatment safety and quality. When the mask has reached its maximum reprocessing limit, discontinue use immediately and dispose of it in accordance with local regulations.
- Regularly clean all components of the mask to prevent bacterial growth and protect your health.
- Do not use aromatic solutions or scented oils (such as eucalyptus or essential oils), bleach, alcohol, or strong-smelling products (e.g., citrus-based cleaners) to clean the mask components.
- After cleaning or disinfection, inspect all parts of the mask for signs of wear, cracks, discoloration, aging, or deformation. Replace any damaged components immediately.
- If the mask has not been used for an extended period, clean and disinfect it before use to prevent dust accumulation and microbial growth, which may affect the safety and quality of therapy.

2 Cleaning and Disinfection

2.1 Cleaning and Disinfection Frequency

		Mask Component	Headgear
Home Environments (Single-patient reuse)	Cleaning Frequency	Daily or after each use	Weekly
Hospital or medical institution Environments (Single-patient/multi-patient reuse)	Cleaning Frequency	Daily/after each use or when switching patients	Weekly/when switching patients
	Disinfection Frequency	Weekly/when switching patients	/

Warning:

- Disinfection must only be performed after cleaning is complete.
- When switching patients, always clean and disinfect the mask first.
- The Headgear does not require disinfection.

2.2 Cleaning (For use in home environments)

This cleaning method is intended for single-patient reuse in a home environment and does not require disinfection.

2.2.1 Mask Component Cleaning

1) Before cleaning, take the mask components apart according to the user manual.
 - N500 series: Detach silicone cushion and mask frame.
 - N300/F300 series: Detach elbow, frame and mask main body.
2) Soak and agitate the mask components (excluding the headgear) in warm water (approximately 30 ℃/86 ℉) using a mild soap for about 10 minutes.
3) Clean the components with a soft bristle brush for about 1 minute, paying particular attention to the small vent holes on the multi-hole elbow.
4) Thoroughly rinse the components under warm flowing drinking quality water (approximately 5 liters

per mask).

5) Inspect each part to ensure that they are visually clean and free of soap residue. If required, repeat washing.
6) Allow the components to dry out of direct sunlight before assembling.

2.2.2 Headgear Cleaning

1) Hand wash the headgear in warm water (approximately 30 ℃/86 ℉) using a mild soap for about 10 minutes.
2) To rinse the headgear, repeatedly squeeze them under warm flowing drinking quality water until it is free of soap residue.
3) Squeeze the components to remove excess water.
4) Rinse well and allow it to air dry out of direct sunlight.

2.3 Cleaning and Disinfection (For use in hospital or medical institution environments)

When reusing this mask among multiple patients in hospitals or medical institutions, the mask components must first be cleaned and then disinfected.

Note: The headgear only requires cleaning and does not need disinfection.

2.3.1 Mask Component Cleaning

Refer to « **2.2.1 Mask Component Cleaning** ».

2.3.2 Headgear Cleaning

Refer to « **2.2.2 Headgear Cleaning** ».

2.3.3 Mask Component Disinfection

Disinfection can be performed using chemical disinfection or thermal disinfection, as follows:

■ **Chemical Disinfection Procedures:**

1) Fully immerse the mask components in a 0.55% ortho-phthalaldehyde solution (e.g., CIDEX OPA) at approximately 20 ℃ (68 ℉) for about 12 minutes.
2) Rinse the mask components with drinking water (approximately 5 liters per mask).
3) Allow the mask components to air dry, out of direct sunlight.

■ **Thermal Disinfection Procedures:**

1) Use a certified hot water disinfection system, soak the disinfectable mask components using a temperature-time combination 1: ISO 15883-1:
 ≥90 ℃ (194 ℉) for ≥5 minutes.
2) On completion, remove the mask components from the hot water disinfection system.
3) Allow the mask components to air dry, out of direct sunlight.

2.4 Maximum Reprocessing Times

Component	Maximum Cleaning Times	Maximum Disinfection Times
Mask Components	180	15
Headgear	24	/

Warning: Exceeding the maximum reprocessing limits may compromise product integrity and safety. Strictly adhere to the specified cycle limits.

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